

REGO X-ray GmbH, Stuttgarter Str. 5 D-86154 Augsburg

Manual Film Processing

rego Processing tank set –3–, all of
PVC for *rego* film hangers up to 35 x 43 cm
consisting of:
2 tanks 20 litres with lid and drain cock
1 tank 30 litres with inlet and overflow
1 drip tray with drain systems
space required: 50 x 56 x 65 cm, net weight 20 kg
Article number _____ 193020

rego Processing tank set –4–, as above but
with additional:
1 tank 20 litres with inlet and overflow
space required: 62 x 56 x 65 cm, net weight 25 kg
Article number _____ 193022

Spare parts and accessories:
Tank 20 litres with lid and drain cock _____ 193024
Tank 30 litres with inlet and overflow _____ 193025
Tank 20 litres with inlet and overflow _____ 193026
Drip tray with drain system for 3 tanks _____ 193027
Drip tray with drain system for 4 tanks _____ 193028
Lid for 20 l tank _____ 193029

Base with steel legs for 3 tanks _____ 193050
Space required: 50 x 48 x 30 cm

Base with steel legs for 4 tanks _____ 193062
Space required: 62 x 48 x 30 cm

Connection set for tank set –3–
consisting of:
2m flexible tube 16 mm inner Ø
2m flexible tube 25 mm inner Ø
4 tube clamps
Article number _____ 193415

Connection set for tank set –4–
consisting of:
3 m flexible tube 16 mm inner Ø
2 m flexible tube 30 mm inner Ø
6 tube clamps, 1 T-piece
Article number _____ 193430

rego Stirring rod, all of PVC
Suitable for *rego* processing tank sets
Stirring rod – grey _____ 196060
Stirring rod – red _____ 196065
Dimensions: 60 cm long, 6 cm wide

Alterations reserved!

