



## X-ray protection goggles



### Fitover Lead Eyewear

Item number 424107

Originally designed to fit over prescription eyewear, this oversized model is also an excellent choice for those requiring good splash protection. The large, clear side shields protect the user from lateral exposure while eliminating the “tunnel” effects of look-alikes.

Material:	Durable Nylon
Lead equivalent:	0,75 mm (+/- 0,05)
Frame Weight:	104 g

Available colour: black (matte)

Available prescription types:  
N/A – for use to wear over existing Rx eyewear

All prints in this brochure are not binding and subject to alterations

Vers.1-dm



**rego** X-RAY GmbH  
**X-RAY ACCESSORIES & MECHANICAL ENGINEERING**  
Stuttgarter Str. 5 86154 Augsburg  
GERMANY  
TEL. +49-821-450 435 0  
FAX. +49-821-450 435 11  
info@rego.de  
www.rego.de

