



**X-ray protection goggles**



### **Panorama Lead Eyewear**

Item number 424108

The wrap frame is constructed of extremely durable Nylon that exhibits superior flexibility while minimizing potential breakage. The molded wrap features provide maximum comfort with a full saddle bridge to distribute weight evenly over the nose while providing excellent peripheral vision with leaded glass side shields.

Frame material:	Durable Nylon
Lead Equivalent:	0,75 mm (+/- 0,05)
Frame weight:	80 g

Colour: tortoise brown

Available prescription types:  
With prescription / non-prescription

Prescription range: + 4,0 and -7,0

All prints in this brochure are not binding and subject to alterations.

Vers.1-dm



**rego** X-RAY GmbH  
**X-RAY ACCESSORIES & MECHANICAL ENGINEERING**  
Stuttgarter Str. 5 86154 Augsburg  
GERMANY  
TEL. +49-821-450 435 0  
FAX. +49-821-450 435 11  
info@rego.de  
www.rego.de

