



X-ray protection goggles



Ultra-light Lead Eyewear

Item number 424106

Manufactured of lightweight Nylon, the Ultralight is designed to offer comfort and provide protection for the entire eye area without the need of an added side shield for lateral exposure. The curved front and temples with rubber tips assures a close secure fit.

Material:	Durable Nylon
Lead equivalent:	0,75 mm (+/- 0,05)
Frame weight:	58 g

Available colours: black or silver

Available prescription types:

Plan:	(non-prescription)
Single:	(low power)
Bifocal:	(with correction)
Prescription range	3,0 to -3,0

All prints in this brochure are not binding and subject to alterations.

Vers. 1-dm



rego X-RAY GmbH
X-RAY ACCESSORIES & MECHANICAL ENGINEERING
Stuttgarter Str. 5 86154 Augsburg
GERMANY
TEL. +49-821-450 435 0
FAX. +49-821-450 435 11
info@rego.de
www.rego.de

