




## X-Ray Film Drying

**Drying Cabinet** for  
20 film hangers up to  
35 x 43 cm (14 x 17")  
for loading from the front.

Equipped with thermostatically  
controlled temperature and a  
maintenance free tangential  
blower.

Metal housing-plastic coated,  
drip tray of stainless steel.

Colour: light-grey



### Technical data:

Supply voltage:	220 V 50 cycles
Current consumption:	600 Watt
Heating:	500 Watt
Timer – adjustable:	0 to 60 min.
Cable connection:	1,5 m Schuco-plug
Capacity:	20 film hangers
Space required:	60 x 50 x 90 cm (h)
Net weight:	25 kg

Article no. 222590